



Celebrating 35 Years of Innovation



CARING SOLUTIONS 2021/2022

Assistive technologies to aid independent living





















3 Melbourne House | Corby Gate Business Park | Corby | Northants | NN17 5JG T +44 (0) 1536 264 869 | F +44 (0) 1536 269719 | E sales@easylinkuk.co.uk

www.medpage-ltd.com | www.easylinkuk.co.uk

(f) http://facebook.com/MedPageLtd @ @MedPageItd

Contents Website ordering 5% discount code MEDCAT21

Visit our websites to view our entire product range

Pages 3 - 5 Page 6 Page 7 Page 8 Page 9	Seizure movement detection alarms Observation cameras (RF & Wi-Fi) Memory prompting & automated messaging GPS location trackers Wi-Fi Connected care at home alarms
Page 10 Page 11 Page 12 Page 13	Medpage long-range transmitter products Medpage MPPL range alarm receivers Medpage MPPL range transmitters Rondish brand alarm receivers
Page 14 Page 15-18 Page 19 Page 20 Pages 21 - 22 Pages 23 - 24	Rondish brand transmitter products Call pendant and special switch products Wireless doorbell solutions - deaf people Gravity change fall detection sensors Door alarms for wander detection Movement sensing alarms (PIR's)
Page 25 Page 26-27	Photoelectric break-beam alarms Floor pressure mat alarms TumbleCare Fall prevention systems
Pages 28 Page 29 Pages 30 - 33	Professional Falls prevention monitor Bed leaving detection alarms
Pages 34 - 36 Page 37	Chair leaving detection alarms SOS emergency help call button dialler

E&OE. We endeavour to maintain 100% accuracy with the content of our publications. We reserve the right to alter pricing and product specification without notice. Please check our websites for any changes.

Trademarks: Medpage, MemRabel and TumbleCare are registered Trademarks owned by Medpage Limited

Orders We accept orders from the public, Local Authorities, NHS and trade resellers (subject to account approval).

© 2021 Medpage Limited T/A Easylink UK. All rights reserved.

EPILEPSY - MEDPAGE SEIZURE MOVEMENT DETECTION ALARMS

Epilepsy Seizure Movement Alarm Medpage MP5V2 High Sensitivity Seizure Monitor



ORDER CODE MP5V2
PRICE £270.00 + VAT
VAT exemption available

- Medpage specialist epilepsy movement detection alarm for Tonic/Clonic seizure movement detection
- Large area patented movement sensor technology for improved movement detection from most bed and seizure movement types
- Tonic/Clonic nocturnal seizure movement detection
- Alarm notifications via radio pager (includes 2 x MPPL pagers)
- Quick and easy to set up
- Suitable for children and adults.
- Suitable for Cots/Cribs with memory foam mattresses, single & double beds

......

The high performance sensor and control software positively identifies seizure movement from a sleeping person, ignoring usual sleeping movements. Verified ongoing seizure movement results in alarm notification transmission to the carer pagers. The sensor (300mm x 200mm x 10mm) has a control to allow adjustment for body weight and bed/cot mattress type. An alarm delay with adjustment from 2-30 seconds allows a person to turn over, or reposition themselves without alarm generation. Included Wi-Fi camera with real time viewing also records seizure activity for viewing on a Smartphone.

Medpage MP5-UT Seizure movement detection alarm Children - Adults

The MP5-UT uses patented sensor technology specifically engineered to detect Tonic/Clonic seizure movement from a sleeping person.

- Detected seizures are notified via the MPPL pagers with tone or vibration alert
- Suitable for babies from age 1 month to adults
- Supplied with 2 carer pagers (1 free) so that both parents can be alerted to detected seizures
- Digital sensitivity control
- Sound sensor detects vocalisations preceding or during a seizure
- Digital alarm delay control virtually eliminates false alarms by ignoring natural sleeping movements such as turning over
- Tamper proof settings & fault reporting to pagers

ORDER CODE MP5-UT PRICE £170.00 (£204.00) VAT Exemption available

MP5-UTB Epilepsy bed movement seizure detection alarm

The MP5-UTB is a reliable nocturnal seizure movement detection alarm for babies through to adults at an affordable price. A portable battery-operated alarm with flashing coloured lights and adjustable audible alarm raises an alert to the carer of detected seizure movement.



- Adjustable sensitivity to suit the needs of the user
- Digital alarm delay control virtually eliminates false alarms by ignoring natural sleeping movements such as turning over
- Digitally stored tamper proof settings
- Inbuilt sound sensor detects vocalisations preceding or during a seizure
- Optional plug in under pillow vibrating pad to assist deaf carers (VIB-PAD)

ORDER CODE MP5-UTB PRICE £170.00 (£204.00) VAT exemption available See page 12 for optional power supply and vibrating pillow shaker for alarm unit

MP2-V2K Medpage seizure movement detection alarm for care homes

The MP2-V2K uses patented sensor technology to detect typical seizure movement from a sleeping person in bed. Carers are alerted to ongoing seizures via a digital message pager (EM300) with tone/vibration alert. Each monitor has a long range transmitter, which is programmed with your chosen message (default for a single alarm system is Medpage Seizure Alarm 1).

- Long range (up to 400M) wireless alarm to carer message pager
- Suitable for use in larger buildings
- · Adjustable sensor sensitivity
- False positive alarm delay control
- Suitable for babies to adults
- Detection of Tonic/Clonic seizure movement
- Vocalisation microphone with alarm to pager
- Can connect to Nurse Call and/or Telecare equipment

Seizure movement monitor including 1 x EM300 pager ORDER CODE MP2V2K PRICE £255.00 (£306.00) See website for more details



EPILEPSY - MEDPAGE SEIZURE MOVEMENT DETECTION ALARMS

MP5-ULTRA Advanced seizure movement detection alarm with detection mode selection and carer alarm pagers

The most advanced epilepsy monitor in the MP5 range. Advanced software algorithms

continually analyse sleeping movements to positively identify seizure activity. Suitable for complex epilepsy and specialist beds including airflow mattresses.



- Detected seizures are notified via MPPL pager (2 supplied) with tone or vibration alert
- Programmed with selectable seizure detection modes
- Digital movement delayed alarm control
- Inbuilt sound sensor to detect vocalisations preceding or during a seizure
- Large area patented monitoring sensor technology offers improved movement detection
- Digital alarm delay control virtually eliminates false alarms by ignoring natural sleeping movements such as turning over
- Tamper proof settings & fail safe alarm to pagers
- Tonic/Clonic and seizure movement with twitching
- Myoclonic cluster seizures
- Airflow mattress mode
- All Medpage seizure movement monitors have 36-hour battery backup
- Suitable for all ages: babies adults
- Nurse call and Telecare systems output connection

ORDER CODE MP5-ULTRA PRICE £439.00 (£526.80) VAT exemption available

All Medpage seizure movement monitors are supplied with a 36 month warranty

Medpage Limited T/A Easylink UK, established 1984, have provided solutions to aid epilepsy care for over 25 years, specialising in the design and manufacture of monitors that detect seizure movements from a sleeping person in bed. To ensure the highest product quality and safety, our company is accredited by BSI to ISO 9001:2015.

At Medpage Limited, we take great care in the design, manufacture and after sales support for our seizure movement alarms. When you purchase an epilepsy seizure movement alarm from Medpage, our highly trained team are readily available for help, advice, and technical support, if needed, to ensure your monitor operates at maximum efficiency. Our technical support is provided at no extra cost.

OBSERVATION CAMERAS (RF & Wi-Fi)

KFM30C Recording Wi-Fi Camera for home and business

View live video monitoring direct to your Smartphone with this simple to install and use home or business WiFi connected camera.

- Audio & video recorded to SD card (sold separately)
- 360° camera rotation via Smartphone APP
- Hi-Def night vision IR LED sensors
- Link with trigger sensors PIR, Door, etc.
- Powered by 5V AC/DC adaptor (included)



Complete camera system ORDER CODE KFM30C PRICE £39.99 (£47.99)

MEDBM-05 Wireless digital baby monitoring camera system

Useful camera system ideal for baby monitoring and patient care.

- Portable rechargeable video monitor
- Battery standby 4-6 hours (typical)
- DC adaptor powered camera
- Wireless range up to 60M (indoors)
- Camera Pan/Tilt from monitor
- High performance night vision



Complete camera system ORDER CODE MEDBM-05 PRICE £83.32 (£98.98)

MEDBM-03 Digital video watch wireless baby monitor

High performance digital wireless video monitor camera with watch style live video viewer.

- VOX vibration alert for camera detected sounds (ideal for deaf carers)
- Wireless range up to 60M indoors (typical)
- Camera night vision viewing
- 2-way voice intercom



Complete camera system ORDER CODE MEDBM-03 PRICE £94.00 (£112.80)


MEMORY PROMPTING & AUTOMATED MESSAGE PLAYBACK

MEMRABEL 2 Digital memory prompting calendar clock for dementia care

This unique dementia clock allows a family member or health visitor to set daily time reminders using voice, picture or video to plan a person's daily routine.

- Create personalised reminders on your phone or PC, or choose from over 100 pre-installed reminders
- File transfer via flash drive or SD card
- Select reminder playback for one-off, daily, weekly
- Choice of time display including traditional analogue
- · Auto seasonal time adjustment
- Powered by mains AC/DC adaptor

ORDER CODE MEMRABEL2 PRICE £99.99 (£119.00)



MemRabel MR3 The smart for solution for routine scheduling

New easy touch screen menus for memory prompt creation with high performance graphics and audio

- Use timed media file playback for daily routine scheduling
- Pre-installed collection of over 150 ready-made reminder prompt files
- Selectable time/calendar display options
- Option to add remote Smartphone APP messaging for £29.99 annually.
- Easy video/slide show player, easy music player
- Visual screen alert for deaf/hard of hearing people
- Touch-screen operation with pass code lock

ORDER CODE: MR3 PRICE £99.99 (£119.00)

Meyro a Arel *2







Membel Remote controlled wireless object locator COL-02

People in the early stages of dementia frequently misplace items around the home. Now, it's very easy to find them. Place a remote beeper on the object (6 included) then simply press the corresponding button on the remote to activate the colour coded beeper.

- Includes 6 beepers powered by CR2032 battery
- Attach to keys, handbags, medicines etc.
- Suitable for the visually impaired 50M location range

ORDER CODE COL-02 PRICE £18.29 (£21.95)



GPS LOCATION TRACKERS - QUICKLY FIND A MISSING PERSON

MMTBV45 Medpage GPS personal location tracker with SOS button, fall sensor and phone calling

- Personal safety for vulnerable and elderly people
- GPS location services via Google maps to Guardian Smartphone
- SMART SOS button for guardian emergency alert notification
- Gravity change sensing fall detection alarm
- Discrete, compact, lightweight, rugged toughened casing
- Waterproof to IP67



MMTBO2 Medpage Micro-GPS tracker with fall sensor & secure Smartphone APP

Rest assured knowing that should a person go missing, they can be located immediately with the Medpage micro tracker.

- GPS location accurate to 5M
- Geo-fence safe zone with breach alarm via APP
- Inbuilt fall sensor sends an alert to emergency contacts should the user experience a fall
- Location by GPS/Wi-Fi hub/GPRS
- Two-way mobile voice calling
- No monitoring or subscription charges
- SOS button calls carer and sends location details
- Includes lanyard and wrist strap

Requires mobile SIM **ORDER CODE MMTBO2 PRICE £58.29** (£69.95) Mobile Roaming SIM offer. Telephone sales for details

MMTBW02 Medpage GPS waterproof tracker watch with heart rate & BP monitor

Lightweight GPS personal tracker watch with fall sensor offers discrete peace of mind to those caring for a loved one who is prone to wandering or falling indoors or outdoors

- GPS tracking accurate to within 5M
- Send personalised reminder alarms to the watch from carer Smartphone - watch vibrates with tone alarm to alert user
- Approximate Heart Rate and Blood Pressure check from carer Smartphone via control APP
- Select preferred time/date display
- SOS function sends location text & 2-way voice calling
- No subscription or monitoring charges
- Waterproof to 1M (for 30 minutes maximum)
- Requires mobile SIM (Optional roaming SIM deal)

ORDER CODE MMTBW02 PRICE £74.99 (£89.99)





Wi-Fi CONNECTED CARE AT HOME ALARMS

RFWFS1K Wi-Fi Smart Plug kit for use with Amazon Alexa, Google Home or with RF pendant transmitter button

The RFWFS1 differs from the majority of Smart Sockets in that it has a built-in radio receiver to allow additional control of the socket with the included pendant transmitter button. The socket can be controlled from a Smartphone, by a Smart Speaker such as Amazon Echo or Google home and by the pendant button, which works independently of a Wi-Fi connection, useful if the internet is not working or outside Wi-Fi range.



Suitable for appliances rated up to 10amp (2200watt)
Control appliances using your Smartphone with Smartlife APP

ORDER CODE RFWFS01K PRICE £29.99 (£35.99)

MEDWFM Wi-Fi connected movement sensor with Smartphone alarm APP

Setting up the MEDWFM takes only a minute or two. Download the user APP to your Smartphone and connect to a home or business Broadband router (2.4GHz only). The user-friendly intuitive APP notifies you if a person walks into the detection beam of the PIR, through the internet to your Smartphone, recording the date and time of the alert. You can check if someone is moving around in the home (parents, elderly relative or potential thief).



- Provides a simple security solution with event date and time record
- Receive an alert if a person enters your home or business premises
- Free to use APP for Android and iOS (no subscription or rental charges)

WiFi movement sensor ORDER CODE MEDWFM PRICE £29.12 (£34.95)

WFSOSB11 Wi-Fi waterproof SOS panic call button for Smartlife APP

- Lightweight waterproof Wi-Fi pendant call alarm for off-site emergency contact
- Wear it with a lanyard or stick it on a wall or hard surface
- APP for Smartphone provides alarm notifications with time/date record (also has 433MHz RF transmitter)
- Quick and easy setup of rechargeable pendant using Tuya Smartlife APP



ORDER CODE WFSOSB11 PRICE: £20.83 (£25.00)

Long-Range POCSAG data message encoded transmitter products

The products detailed on this page are featured in various kits in the following pages. An unlimited number of transmitters can be used with single or multiple pagers. The messages are factory programmed so you need to provide us with your messages when placing an order for additional or multiple products. Single items are shipped with the appropriate default message. Some items are not featured in this brochure.



POC-840DC Quick and easy to install door & window alarm with on/off switch. Programmable with message up to 40 characters.

ORDER CODE POC-840DC PRICE £62.00 (£74.40)

400M range Power 1 x CR2032 Low battery warning



POC-940 Waterproof (IP67) call pendant with wrist strap and lanyard included. Programmable with message up to 40 characters.

ORDER CODE POC-940 PRICE £58.00 (£69.60)

400M range Power 1 x CR2450 Transmit confirmation LED - Low batt alert



NMDTX2 Multi-function message programmable POCSAG transmitter with configurable N/O N/C input. ORDER CODE NMDTX2 PRICE £56.00 (£67.20)

400M range Power 2 x AAA Bed/chair & floor pressure mat modes.



EM300 Digital data message pager with tone/vibrate alert. Stores last 20 messages with time & date. **ORDER CODE EM300 PRICE £70.00 (£84.00)**

Sturdy belt clip Power 1 x AAA Time/date calendar Time alarms



EM400 Digital data message pager with powerful vibration and tone alert. Top view vertical display. **ORDER CODE EM400 PRICE £61.00 (£73.20)**

400M range Power 1 x AA LCD backlit display Message review



346BC Long range call button with inertia and gravity fall sensor. IP67 waterproof. **ORDER CODE BC346 PRICE £78.00 (£93.60)**

400M range Power 1 x CR2450 Message programming Transmit & fall LED



EM1600U Signal repeater for all Medpage long-range transmitter products (this page). Extend range up to 500M. Suitable for Care Home and security applications ORDER CODE EM1600U PRICE £147.00 (£176.40)

Range up to 500M Power 12V DC @ 1A



MP-900 Radio frequency to Wi-Fi adaptor. Receive alerts to pagers and Smartphone Via free APP ORDER CODE MP-900 PRICE £115.00 (£138.00)

Requires home or business Broadband. Perfect combination for security and home care.

MEDPAGE WIRELESS ALARM RECEIVERS & ACCESSORIES

Medpage manufacture the alarm receivers detailed below to offer choice for alerting carers. Some people prefer a pager, others a tabletop alarm. The CTMV allows recording of personal alarms or reminders, triggered by remote transmitter. Where you see the Medpage logo in the following pages, it indicates the products are compatible with these alarm receivers.



Portable alarm pager with belt clip. Robust in design and very easy to use. Compatible with all MPPL transmitters.

Power: 2 x AA Tone/vibration alarm Adjustable volume For up to 9 devices Typical 100M range



Portable 5-Channel alarm station with tone & individually coloured flashing light

alarms. Compatible with all MPPL

ORDER CODE MPPL2 PRICE £31.00 (£37.20)

transmitters.

Power: 3 x AA Tone/flashing alarm Adjustable volume For up to 5 devices Typical 100M range

ORDER CODE NMDRX £37.49 (49.99)



Plug-in vibrating under-pillow shaker for use with the NMDRX alarm station for the deaf and hard of hearing.

Rated at 300mA Cord length 2M 3 Volt (MAX 5V) 3.5mm Jack plug





Portable voice alarm receiver with 3 recordable voice channel alarms of up to 15 seconds. Relay output allows any Medpage transmitter to connect to Nurse Call or Telecare alarm.

Power: 3 x AAA Voice recorder alarm Adjustable volume For up to 3 devices Typical 100M range

ORDER CODE CTMV PRICE £35.83 (£42.99)



Medpage signal range extender compatible with all MPPL/Medpage transmitters.

ORDER CODE NMDTX2 PRICE £56.00 (£67.20)

UK Plug top AC/DC power adaptorPlease enquire for correct product model

Universal operating mode transmitter with signal range extender for Medpage standard and long range transmitters.

MEDPAGE BRAND WIRELESS TRANSMITTER ALARMS

All products featured in this section are additional transmitters for use with the MPPL pager, NMDRX 5-Channel alarm or CTMV 3-Channel alarm. They can be ordered separately where extra or mixed sensors are required. All products have an approximate transmitter range of 100M on 433MHz. Use the product code to search our website.

Universal transmitter with N/O N/C volt free trigger input ORDER CODE NMDTXM PRICE £33.32 (£39.98)





Waterproof push button call pendant with lanyard and wrist strap ORDER CODE MPPL-ERTX Price £14.00 (£16.80)





Wireless door/window alarm transmitter (MPPL pager only) ORDER CODE MED-DCT PRICE £11.66 (£13.99)





Wireless doorbell button transmitter ORDER CODE NMD-DB1 PRICE £8.28 (£9.94)





Sound activated transmitter for baby monitoring and voice help call ORDER CODE MPPL-SAT Price £21.00 (£25.20)





Battery powered PIR with alarm pause button and selectable audio or transmitter alarm. Signals all Medpage alarm receivers





ORDER CODE TUMPIR21 Price £26.00 (£31.20)

Wrist worn fall sensor with touch sensor call button ORDER CODE MPPL-FA6 PRICE £27.45 (£32.94)





Washable cotton enuresis sensor with alarm transmitter - Long range version available. Please enquire for details. ORDER CODE NMDTXM-ES01 PRICE £61.72 (£74.06)



PIR Movement sensor transmitter with external ON/OFF switch ORDER CODE MPPL-MS PRICE £28.00 (£33.60)





RONDISH BRAND WIRELESS ALARM RECEIVERS

Medpage is the UK distributor for Rondish nursing care products. The products featured on this page are alarm receivers, which can be used individually or in combinations for individual or multiple room/ward monitoring care systems. Where you see the Rondish logo in the following pages, it indicates the products are compatible with these alarm receivers.

ILB-11 Wall mounting audible & visual alarm



- Can pair with up to 32 transmitter products
- Audible alarm with volume adjustment
- Bright flashing alarm LED's
- Integral signal booster for all paired transmitters

ORDER CODE ILB11 PRICE £112.00 (£134.40)

CMD-11c Desktop alarm station with signal booster



- Desktop wireless caller display alarm (99 devices)
- Integral signal booster for paired devices
- Mains powered with battery backup (10 hours)
- Multiple call display call review

ORDER CODE CMD-11 PRICE £105.00 (£126.00)

PAG-11C Caller display pager with audible & vibration alert alarm



- Tone/vibration alert carer pager with caller display
- Can pair with up to 99 Rondish transmitters directly
- Sensor type and caller ID assigned during setup
- Powered by 2 x AA batteries (typical 6 months)

ORDER CODE PAG-11C PRICE £31.00 (£37.20)

CMU-02 Wireless nurse call alarm station for up to 999 devices



- Mains powered with battery backup (4 x D-Cell)
- Call and system fault reporting
- Call silence button with delayed alarm repeat
- Displays call type, zone and room/ward ID
- Fully monitored device connection status

ORDER CODE CMU-02 PRICE £205.00 (£246.00)

RONDISH BRAND ADDITIONAL/REPLACEMENT ACCESSORIES

BMS-03 Non-radio alarm controller for SEN-1/SEN-2/FMAT-2 alarm pads ORDER CODE BMS-03V2 PRICE £38.00 (£45.60)



BMDV-11 Non-radio alarm controller for SEN-1/SEN-2/FMAT-2 alarm pads Magnetic pull-away clothing clip alarm and nurse call connection ORDER CODE BMDV-11 PRICE £54.12 (£64.94)



DMS-02 Replacement/additional door sensor for DoorWatcher egress alarm ORDER CODE DMS-02 PRICE £248.00 (£297.60)



TXP-02R Replacement/additional resident wristband for DoorWatcher alarm **ORDER CODE TXP-02R PRICE £46.00 (£55.20)**



TM11(V3) Replacement transmitter for cordless bed/chair pads ORDER CODE TM11(V3) PRICE £31.45 (£37.74)



WCP11 Waterproof pull-cord transmitter - use with any Rondish receiver **ORDER CODE WCP-11 PRICE £56.00 (£67.20)**



DRB-11 Wireless call reset for DoorWatcher system **ORDER CODE DRB-11 PRICE £68.00 (£81.60)**



SEN-1 Bed sensor for BMS-03, BMDV11, BTX-03M **ORDER CODE SEN-1 PRICE £28.00 (£33.60)**



SEN-2 Chair sensor for BMS-03, BMDV11, BTX-03M **ORDER CODE SEN-2 PRICE £24.00 (£28.80)**



Wire-Free bed and chair occupancy sensor pad with integral transmitter Bed pad ORDER CODE CBP-01 PRICE £70.10 (£84.84) Chair pad ORDER CODE CCP-01 PRICE £65.80 (78.96)

New dual sensor bed edge detection pressure pad with integral wireless alarm transmitter. Compatible with all Rondish brand alarm receivers. Designed for use with the MPCSA11 falls system.

ORDER CODE WPAS10 Price £61.53 (£73.84)

CALL PENDANT AND SPECIAL SWITCH PRODUCTS

MPPL-ERTXSET Waterproof call pendant with MPPL2 pager

Easy to press pendant alarm allows the user to raise an alert to a carer's portable pager, simply by pressing the pendant button.

- Simple wireless SOS call bell solution
- · Ideal for the elderly or infirm
- Tone and vibration alert pager
- Pager can pair with up to 9 transmitters
- Pendant is supplied with lanyard and wrist strap
- · Portable, lightweight design
- All batteries included



NMDRX-ERTX Waterproof wireless call pendant with alarm station

The call pendant button signals the alarm station, which is equipped with coloured bright flashing LED's. Useful for hard of hearing or deaf carers.

- Simple to use and set up call solution
- 5-Channel independently coloured flashing LED's
- Output for under pillow vibrating pad
- Front facing and ceiling projected flashing LED's
- Long battery life with low battery warning
- Audible alarm tone with 80-90 dB selection
- Choose your preferred alarm light colour per channel



Call button with alarm ORDER CODE NMDRXERTX PRICE £51.49 (£61.79)

PAG31ERTXKIT Waterproof call pendant with 3-channel alarm player

Press the waterproof call button to alert a carer via the pager that help is needed.

- Wireless no-fuss call button with alarm pager
- Ideal for shower and wet room help call to carer
- 3-Channel pager can monitor 3 transmitter devices
- Pager alert notification by audible and vibration alarm
- Pager powered by 2 x AA batteries
- Pendant transmitter range 50-80M



Pendant with pager ORDER CODE PAG31ERTXKIT PRICE £42.20 (£50.64)

CALL PENDANT AND SPECIAL SWITCH PRODUCTS

MPPL-PSKIT Soft touch low tactile pillow switch with MPPL pager

Face/head activated switch attaches to pillow to enable a person with limited mobility to activate a help call to a carer pager.

- Soft and comfortable washable switch cover
- Attaches to pillow with included clips
- Signals carer tone/vibration alarm pager
- Transmitter range up to 100M
- Suitable for stroke victims & disabled people
- All batteries included



Pillow switch kit ORDER CODE MPPL-PSKIT PRICE £164.00 (£196.80)

CTMV-PSKIT Soft touch low tactile pillow switch with voice alarm

Face/head activated switch attaches to pillow to enable a person with limited mobility to activate a help call to a carer voice player.

- Use personalised voice recordings as an alarm
- Record up to 3 separate call messages (additional transmitters required)
- Suitable for child carers or people with mild dementia
- Record an instruction for playback when the button is pressed
- Suitable for stroke victims & disabled people
- All batteries included

Pillow switch kit ORDER CODE CTMV-PSKIT PRICE £169.32 (£203.18)

NMDRX-PSKIT Soft touch low tactile pillow switch with alarm station

Face/head activated switch attaches to pillow to enable a person with limited mobility to activate a help call to a carer.

- 5-Channel alarm receiver with 5 selectable alarm colours
- Suitable for hard of hearing/deaf carers
- Audible/flashing light alarm station
- 30 Second or looped audible/flashing light alarm
- Transmitter range up to 100M
- Suitable for stroke victims & disabled people
- All batteries included

Med Page

Pillow switch kit ORDER CODE NMDRX-PSKIT PRICE £170.99 (£205.19)

CALL PENDANT AND SPECIAL SWITCH PRODUCTS

MPPL-JSK Large surface area call button switch with carer alarm pager

For people that find it difficult or impossible to operate small call pendant buttons, the jelly switch provides a perfect solution. Available with a choice of wireless carer alarm responders.

- Large surface button can be activated with a minimal amount of pressure
- Suitable for patients with poor dexterity
- Image can be placed under the cover to assist patients with dementia or learning difficulties
- Wireless signal to tone/vibration alert alarm pager
- All batteries included

Jelly switch kit with pager ORDER CODE MPPL-JSK PRICE £82.00 (£98.40) Switch with transmitter ORDER CODE MPPL-JS PRICE £51.49 (£61.79)

PAG31-JSK Large surface area call button switch with pager alarm

This jelly large button area switch kit includes a 3-Channel caregiver pager with a universal transmitter and Jelly Switch.

- Easy to press large call button suitable for people with Arthritis, Gout, Stroke or physical disability
- Caregiver pager with audible and vibration alert.
- Transmission range 100M
- Quick and easy call button solution



Med M

Jelly switch kit with pager alarm ORDER CODE PAG31-JSK PRICE £78.71 (£94.44) Switch with transmitter ORDER CODE MPPL-JS PRICE £51.49 (£61.79)

NMDRX-JSK Large surface area call button switch with carer alarm station

This jelly switch kit uses the NMDRX call station with audible/flashing light alarm.

- 5-Channel alarm receiver with 5 selectable alarm colours
- Suitable for hard of hearing/deaf carers
- 30 second or looped audible/flashing light alarm
- Transmitter range up to 100M
- Suitable for debilitation conditions: Stroke, MS, Parkinson's etc.

Jelly switch kit with alarm station ORDER CODE NMDRX-JSK PRICE £88.49 (£106.18) Jelly switch with transmitter ORDER CODE MPPL-JS PRICE £51.49 (£61.79)



CALL PENDANT FOR SINGLE AND MULTIPLE USERS

PAG-11-ULC11 Waterproof call pendant button with caller display alarm pager

The PAG-11 pager is equipped with tone and vibration alert alarms. Transmitters are assigned to a user ID (up to 9), for pager caller ID display

- Suitable as a single or multi-user call pendant alarm
- Waterproof transmitter with soft touch call button
- Call review button alarm cancel button setup menu
- Low battery warning for call button and pager
- PAG-11 is compatible with all Rondish transmitters



Call pendant with pager ORDER CODE PAG-11-ULC11 PRICE £54.95 (65.94)
Call pendant only ORDER CODE MPPL-ERTX PRICE £14.00 (£16.80)

POC-940K Long range waterproof call pendant with digital message pager

This call pendant has a typical transmission range of up to 400M. It is compact, lightweight and can be used with either a lanyard or wrist strap, both supplied. The data message pager has tone and vibration alert.



- Suitable for domestic and professional use
- Unlimited call pendants can be monitored by a pager
- . Long battery life (up to 6 months) with low battery warning
- Suitable for larger and stone built buildings
- Stores call alerts with date and time

Transmitter and pager kit ORDER CODE POC-940K PRICE £128.00 (£153.60)
Pendant transmitter only ORDER CODE POC-940 PRICE £58.00 (£69.60)

RONWC2 Waterproof pull cord alarm for disabled toilet with emergency call alarm

This is a complete disabled toilet alarm system for use in public and commercial buildings. The battery powered emergency call alarm is wall mounted and has an audible and flashing light alert when signalled by pulling the cord.

- Wireless pull cord transmitter (Range 60M O/F)
- 3M cord length with red pull triangles
- On-device SOS call & wireless reset button
- Install for CQC and DDA compliance



Rondish

Complete system ORDER CODE RONWC2 PRICE £168.00 (£201.60)
Additional pull cord ORDER CODE WCP-11 PRICE £56.00 (£67.20)

WIRELESS DOORBELL BUTTON WITH ALARM CHOICES

NMDRXDB Wireless doorbell button with flashing light alarm receiver

Simple to install and use flashing light doorbell system, ideal for the deaf and hard of hearing.

- Battery operated portable tabletop alarm
- 5-Channel with selectable flashing alarm colour
- Receiver can pair with up to 5 transmitters
- Output for optional under pillow vibrating pad
- 30 second or looped alarm
- Secure doorbell fixings supplied
- All batteries included



Flashing light doorbell kit ORDER CODE NMDRXDB PRICE £37.49 (£44.99)
Doorbell button only ORDER CODE NMD-DB1 PRICE £8.28 (£9.94)

MPPL-DBK Wireless doorbell button with tone/vibration alarm pager

Simple to install and use radio paging doorbell system, ideal for the deaf and hard of hearing.

- Doorbell transmitter button activates pager alarm
- Transmission range up to 80M
- Pager can pair with up to 9 transmitters
- Sturdy pager belt clip
- You can add an NMDRX flashing light (to identify active doorbell signal) alarm or additional pagers





Doorbell with pager ORDER CODE MPPL-DBK PRICE £31.00 (£37.20)
Doorbell button only ORDER CODE NMD-DB1 PRICE £8.28 (£9.94)

PAG31DB Wireless doorbell button with alarm pager

Our best priced wireless doorbell solution for deaf people. Pager vibrates when doorbell button is pressed. You can also add visual alarm NMDRX. Quick and easy to install. Ideal for home use, students in halls or at work etc.



- 3-Channel pager can be used with 3 doorbell buttons
- Pager 2 x AA batteries (included)
- The doorbell button in this kit is free

Doorbell with Vibration alert pager ORDER CODE PAG31DB PRICE £28.20 (£33.84)
Doorbell button only ORDER CODE NMD-DB1 PRICE £8.28 (£9.94)

GRAVITY CHANGE FALL DETECTION SENSORS WITH ALARM

BC346-EM300 Waterproof fall detection bracelet with long range call button

This kit includes a high specification inertia/gravity fall sensor (software adjustable) with convenient wrist strap and call button and offers at home and away from home (with optional MP-900 iOT internet hub) fall notification and includes a high performance tone/vibrate radio pager. With a SOS and fall detection transmitter range of up to 400M the kit is suitable for domestic and Care Home use. The optional Wi-Fi connected control hub provides off-site alarm alerts and setup of any sensors added to the kit including: door alarms, bed & chair occupancy sensors. Kit includes: 1 x fall sensor bracelet, 1 x EM300 pager.

ORDER CODE BC346-EM300 PRICE £139.00 (£166.80)
ORDER CODE MP-900 PRICE £115.00 (£138.00) (optional)



MP900

TUM31TXP11 Gravity activated fall detection alarm with pager

For people caring for a person prone to falls at home or in care this kit will provide a reliable alert to a fall. Additionally, the call button can be pressed to signal for help. Ideal for use when showering etc.

- Convenient SOS button with integral fall detection sensors
- Provides rapid alert to a carer of a potential fall
- Rugged alarm pager with tone/vibration alert
- Flashing active alarm channel LED
- Pager suitable for deaf and hard of hearing carers
- Pendant waterproof to IP67

ORDER CODE TUM31TXP11 PRICE £59.88 (£71.86)
Fall alarm only ORDER CODE TXP-11 PRICE £32.28) (£30.94)

NMDRX-FA6P Wrist worn gravity activated fall detection alarm

This kit uses the NMDRX call station with audible/flashing light alarm. You could add the CTM3-CT3 bed alarm to provide daytime fall and night time out of bed alert.

- 5-Channel alarm receiver with 5 selectable alarm colours (example: fall alarm RED bed alarm BLUE)
- · Suitable for hard of hearing/deaf carers
- 30 second or looped audible/flashing light alarm
- NMDRX power: 3 x AA batteries or optional AC adaptor

Or Page

Fall alarm with NMDRX kit ORDER CODE NMDRX-FA6P PRICE £64.94 (£77.93) Fall alarm only ORDER CODE MPPL-FA6 PRICE £27.45 (£32.95)

WIRELESS DOOR ALARMS FOR WANDER DETECTION - SECURITY

TUMPAG31DCT Medpage wireless door opening alarm with carer alarm receiver

Simple to install door security solution. When a monitored door or window is opened, an alarm signal is transmitted to the alarm pager. The pager can monitor 3 separate doors with device ID shown by channel LED

- Dementia and elderly wander detection
- General door and window security
- Door alarm call/panic button
- Sound and visual alert for activated sensor transmitters
- 3-Channel tone/vibration alert pager
- Typical range 100M (line of sight)
- Pager powered by 2 x AA batteries DCT CR2032 battery

Door transmitter with pager ORDER CODE TUMPAG31DCT PRICE £39.26 (£47.11)

MPPL-DCT Medpage wireless door opening alarm with carer alarm pager

Where a mobile alert to an opening door is required, this kit provides a perfect solution. When the door/window is opened an alarm signal is transmitted to the alarm pager.

- Dementia and elderly wander detection
- General door and window security
- Alarm pager powered by 2 x AA batteries
- Tone/vibration alert with volume selection
- Door alarm powered by 1 x CR2032 battery
- Typical transmitter range 60M



Door transmitter with pager ORDER CODE MPPL-DCT PRICE £42.66 (£51.19)

Emergency call pendant with qsm telephone auto dialler ELK4-PO

- Mobile network emergency phone dialler
- Set up of dialler by text messages
- Receive text and phone call alarm notification
- Includes free wireless door sensor (while stocks last)
- Waterproof call pendant with lanyard
- Two-Way voice intercom & on device SOS call
- Monitors home temperature with Hi/Low alarms
- Requires mobile SIM (see 1P Mobile offer)

ORDER CODE ELK4-PO PRICE £40.00 (£48.00)



WIRELESS DOOR ALARMS FOR WANDER DETECTION - SECURITY

PAG-11C-DCT Wireless door alarm with caller display tone/vibrate pager

This kit provides a simple door alarm solution for single or multiple doors. The caller display will show a user ID to indicate door location.

- Suitable for domestic, Care Home and Hospital use
- Use to detect wandering people exiting internal or external doors
- Can integrate with all Rondish wireless care aids
- Pager powered by 2 x AA batteries
- Up to 9 door alarms (or combined devices) can be paired with the pager
- Typical transmitter range 60M. See signal boosters.

Door kit ORDER CODE PAG-11C-DCT PRICE £42.66 (£51.19)

MEDWFD Wi-Fi connected door alarm with monitoring APP

The MEDWFD is a useful door/window alarm sensor, which connects through your home or business Broadband (2.4GHz only) and is configured via a Smartphone APP (Android or iOS). The APP provides a comprehensive selection of features, including choice of alarm notifications and sensor labelling.

- APP records date & time of door events
- · Can also be used on food cupboards and fridges
- Beneficial in the care of people with dementia or cognitive impairment for wander detection
- Can be used with multiple sensors for home and business protection

Wi-Fi door/window sensor ORDER CODE MEDWFD PRICE £24.95 (£29.94)

POC-DCK Long range (400M) door alarm transmitter with digital message pager

When a door or window alarm transmitter is required with longer transmission, this kit provides a good solution. The compact transmitter (POC-840DC) is programmed to send location messages to the EM300 message displaying pager. Unlimited transmitters can be used with a single or multiple pagers.

- Use for wander detection & general door/window security with door location message
- Suitable for single or multiple alarm transmitter systems
- Suitable for domestic, Care Home and Hospital use
- Factory programmed default message 'Door Alarm 1'
- Tone/vibrate alert pager powered by 1 x AAA battery

Long range door alarm kit ORDER CODE POC-DCK PRICE £132.00 (£158.40)



MOVEMENT SENSING ALARMS - BED LEAVING - WANDERING

MPPL-MSKIT PIR movement sensor transmitter with pager

This easy to install movement alarm kit can be used for general security. It is also beneficial to aid care for the elderly and those living with dementia. When the sensor detects movement an alarm signal is transmitted to the pager.

- Use for bed leaving detection and wander alerts protect stairwells and doorways
- Long battery life for PIR & pager (6-12 months)
- Transmitter range up to 100M
- External power on/off switch for PIR
- Tone/vibration alert pager



DRXTMSKIT Movement detection PIR wander detection kit with portable alarm

Detect an elderly patient or a person with dementia wandering during day or night with this easy to use movement detection alarm kit.

- An alert is transmitted to a carer when a person walks through a doorway, enters a room, exits the home or enters a high-risk area
- Includes NMDRXT 5-Channel carer alarm receiver
- Quick and easy installation
- Narrow PIR detection angle allows use as a break beam alarm
- Quick and easy to set up use straight out the box

PIR with NMDRXT receiver ORDER CODE DRXTMSKIT PRICE £63.00 (£75.60)

TUMMSKIT Narrow beam PIR movement sensor with alarm pager



Motion detection sensor with alarm pager kit. The motion sensor has a switch to select audible alarm or integral 433MHz transmitter. On movement detection the selected alarm type is activated, pager or audible tone.


- Movement sensor with pager alarm
- Perfect to warn a carer of elderly wandering or bed exit
- Narrow beam PIR sensor
- 3-Channel pager with tone/vibration alert
- Cargiver movement alarm pause button

PIR with Pager ORDER CODE TUMMSKIT PRICE £54.00 (£64.80)

LONG RANGE (400M+) POCSAG MESSAGING SOLUTIONS

POC-PIR Message programmable Movement Sensor

When multiple Movement Sensors are required in larger buildings such as Care Homes, Hospitals or buildings known to have signalling difficulties, the POC-PIR is the perfect solution. Furthermore, messages can be programmed to enable location, device, or zone ID to be read from our digital message display pagers (EM300/EM400). Ideal for wander detection, bed leaving detection and general security use.



- Long range transmission 400M+
- Radio frequency: 433MHz
- Low standby power consumption
- Side power on/off switch
- Dimensions: 107 x 59 x 39mm
- Weight: 40g
- Power: 2 x AA batteries
- Supplied with wall mount bracket
- Status LED for power/transmission/signal
- Compatible with EM300 & EM400 pagers
- Compatible with MP900 internet hub
- Supplied with standard message "PIR ALARM" and configured to standard pager radio address code.

Programming software download available for **Windows**



ORDER CODE POC-PIR PRICE £46.40 (55.68)

MULTI-FUNCTION LONG RANGE UNIVERSAL TRANSMITTER 400M+

The NMDTX2 is designed to operate with all Medpage pressure mat sensors and can be programmed to transmit messages for display on the EM300/EM400 message display pager.

Selectable operating modes

- Bed occupancy detection
- Chair occupancy detection
- Floor pressure mat
- MPPL signal repeater
- Alarm input trigger socket N/O N/C
- Transmitter range 400M+
- Message programmable
- Frequency 433MHz
- Can connect to MP900 iOT hub for Smartphone notifications
- Status LED power/transmit
- Power 2 x AAA/AC adaptor



ORDER CODE NMDTX2 PRICE £56.00 (67.20)

MOVEMENT SENSING BREAK BEAM ALARMS

Solar powered break-beam photoelectric intruder alarm MED-PEB



- Solar powered photoelectric infrared beam alarm
- Built in radio transmitter to signal Medpage alarm receivers (MPPL, TUMPAG31, NMDRX
- Quick and easy beam alignment up to 60M
- No need for external power sources
- Alternative to PIR movement sensors
- Alarm generated only on direct break-beam
- Invisible infrared dual beams
- No false alarm from small rodents or birds
- Easy beam alignment using fine tuning controls

ORDER CODE MED-PEB PRICE £150.25 (£180.30) PRICE FOR BEAM SET ONLY

This beam set provides a quick and easy security solution indoors or outdoors (IP67). An intruder or larger animal passing through the beam will generate an alarm transmission with a range up to and beyond 100 Metres. The Medpage alarm receivers MPPL, TUMPAG31, NMDRX provide wire-free notification of a security breech.

FT47C-DTXK Photoelectric break beam sensor alarm kit with long range transmitter and digital message pager

The FT47C-DTXK is a compact photoelectric beam set suitable for indoor and outdoor use (shower proof only) to create an invisible narrow beam security solution. An object, person or animal walking through and breaking the beam results in an alarm transmission to the EM300 data message pager.

- Long range alarm transmission (up to 400M)
- Beam range (strong sunlight) 8-12M (low sunlight) 15-17M
- Programmable alert messages (when used with data pager EM300)
- Working voltage IR transmitter 2 x AA Alkaline batteries
- Typical battery life 12 months 10 ops per day
- Working voltage IR receiver: 12V DC power adaptor or battery
- IP rating: IP55 outdoor use (shower proof)



ORDER CODE FT47C-DTXK

PRICE £151.00 (£181.20)

FLOOR PRESSURE MAT ALARMS

MPPL-PMK Floor pressure mat alarm with radio pager

Our best-selling floor pressure mat alarm is ideally suited to the domestic user. The mat is easily concealed under a rug or carpet next to an exit door or the user's bed. When the user steps on the mat, the carer is notified via the MPPL pager.

- 100M transmission range
- 3 alert options: vibrate only, vibrate and low tone or vibrate and high tone
- Robust, lightweight pager with a sturdy belt clip
- Low battery warning for uninterrupted user monitoring
- Pressure mat measures 750mm x 450mm with a 3M connection lead

Complete kit ORDER CODE MPPL-PMK PRICE £75.00 (£90.00)

TUMFMMPPLK TumbleCare large heavy duty non-slip floor pressure mat alarm kit

The non-slip heavy duty pressure mat is positioned on the floor next to the patient's bed. Standing on the mat results in alarm transmission to the caregiver pager. Suitable for home care, Care Home and Hospital use. Can also be used at entrance doors and for general security.

- Mat measures 71.44cm x 60.96cm
- Dual surface carpet grip and rubberised surface
- Includes nurse call transmitter (BTXMED04)
- Includes caregiver pager (TUMPAG31)
- Suitable for single or multi-patient monitoring (up to 3)
- 3-Metre sensor connector cable

Complete kit ORDER CODE TUMFMMPPLK PRICE £129.81 (£155.77)

PM-POCKIT Floor pressure mat alarm with long range (400M) message transmitter & data message pager

The PM-POCKIT is recommended for larger buildings, stone-built properties or buildings with known signalling difficulties. The pressure mat (750mm x 450mm) can be positioned under rugs or carpets and will transmit an alarm signal when a person stands on it.

The kit comprises of NMDTX2 long range POCSAG transmitter, 1 x CT-3 pressure mat and 1 x EM300 pager. An unlimited number of long range devices can be configured into a multiple care system. Single kit default message "Pressure Mat Alarm".

Complete single kit ORDER CODE PM-POCKIT PRICE £142.17 (£170.60)





FLOOR PRESSURE MAT ALARMS

DTXTUMO4K TumbleCare Heavy duty non-slip floor pressure mat with long range transmitter (400M) & digital data message pager

This kit is recommended for larger buildings, stone-built properties or buildings with known signalling difficulties. The pressure mat (91cm x 61cm) can be positioned on top of rugs, carpets, or directly on hard floors and will transmit an alarm signal when a person stands on it.

The kit comprises of NMDTX long range POCSAG transmitter, 1 x MED-04FP pressure mat and 1 x EM300 pager. An unlimited number of long range devices can be configured into a multiple falls care system. Single kit default message "Pressure Mat Alarm".

The state of the s

Single system kit ORDER CODE DTXTUM04K PRICE £182.00 (£218.40)

NMDRX-CTMFP Floor pressure mat alarm with transmitter & portable alarm

This kit includes the CT-3 pressure sensor, NMDRX audible/flashing light alarm and the CTM-3 voice alarm transmitter. A warning message can be recorded for playback when a person stands on the mat. A signal is transmitted to activate the NMDRX.

- Dual purpose alarm voice prompt and carer alert
- Use for dementia care wander detection.
- NMDRX can pair with up to 5 sensors to create expanded care systems



Complete kit ORDER CODE NMDRX-CTMFP PRICE £102.99 (£123.59)

WT-ELB Wi-Fi connected floor pressure mat alarm with Smartphone alerts

For families concerned about elderly relatives living alone, the WT-ELB floor pressure mat alarm will provide peace of mind and valuable time/date stamped information of daily living habits. The powerful free Smartphone APP allows you to set times when the sensor is active and choose the type of alert you want to receive to your phone. When activated, the Wi-Fi controlled sensor will send an alarm signal to your Smartphone via your home/business Wi-Fi connection.

- Wander detection fall prevention
- Controller powered by 2 x AAA batteries
- Under carpet/rug pressure sensor (750mm x 450mm)
- Pair with home/business broadband (2.4GHz only)

Complete kit ORDER CODE WT-ELB PRICE £86.66 (£103.99)

Tumble Care FALL PREVENTION SYSTEMS FOR PROFESSIONALS

BTXMED04 pressure sensor alarm transmitter (433MHz)

- Pressure sensor pad alarm controller
- Wireless alarm signal to carer receiver
- Tough durable design
- Supplied with rubber protection jacket
- Selectable alarm response (instant, delayed 15 minutes, delayed 30 minutes)
- Suitable for individual or multi-patient systems
- Nurse call peer push input socket
- Central illuminated patient call button
- Alarm reset button
- Sensor pad strain relief protector

ORDER CODE BTXMED04 PRICE £41.21 (£49.45)

MED-04FP Large heavy duty non-slip floor mat

MED-04FP Floor pressure sensor mat supplied with connector adaptor to enable compatibility with all TumbleCare, Medpage/Easylink monitor transmitter products.

Dimensions: 91.44cm x 60.96cm

Virtually indistructible

Cable length: 3M

Connector type: 6.3mm mono jack plug

Chlorine resistant to 10pp1000



ORDER CODE MED-04FP PRICE £61.00 (£73.20)

TUMPAG31 New - TumbleCare 3-Channel caregiver alarm notification pager



The new TUMPAG31 3-channel pager by Medpage provides reliable carer notifications with audible tone & vibration alert. Each channel has an active alarm indicator LED. Monitor up to 3 bed or chair alarms, door or movement sensors.

ORDER CODE TUMPAG31 PRICE £28.20 (£33.12)

TumbleCare is a new brand of fall detection/prevention products by Medpage, launched 2020/2021. With over 30 years worth of experience within the Healthcare sector, our skills and knowledge are focused into our new range of products designed to deliver affordable, reliable performance for home and professional care providers.

MPCSA11 INTELLIGENT SENSOR ALARM CONTROL TRANSMITTER

Bed and chair falls wireless sensor pad alarm controller

The MPCSA11is a radio linked sensor monitor transmitter, designed to provide the essential features required for a modern patient falls monitoring system while maintaining ease of use and setup. The CSA11 is compatible with cordless and cabled bed, chair and floor pressure mat sensors.

A feature unique to the CSA11 is last sensor detection alarm reset. This feature can help to track a person leaving their bed and moving to a chair without an alarm generation. Additionally, a movement sensor can be added to create a comprehensive mobility monitoring system with a last detection timeout alarm to a carer radio pager or telephone emergency dialler to alert off-site carers to a potential fall.



- Suitable for single or multiple patient bed/chair exit with pager alarm notification
- monitor 3 cordless sensors with RF link alarm output to nurse call system
- Peer push switch and nurse call connection socket
- Turn patient alarm pause button
- Monitoring for up to 3 wireless sensors: bed, chair, PIR, door alarm

ORDER CODE MPCSA11 PRICE £78.10 (93.72) FULL SPECIFICATIONS

Competitively priced kits available for professional care. Contact sales for details

New Bed exit sensor with dual pressure sensors

Dual sensor bed edge detection pressure pad for use with the MPCSA11 falls system. Edge sensors detect a person sitting or rolling to the edge of the bed. The central sensor detects body weight from the middle of the bed.

- Suitable for single & double size beds
- Supplied with clips to secure sensor to mattress
- Central and bed edge body weight detection
- Chlorine resistant sensor material
- System transmitter concealed in mat pocket
- Mat dimensions L 82cm W 26cm D 5mm
- Bed exit alarm transmission notification to alarm receiver (PAG-11, TUMPAG31, CMD11, MPCSA11)
- Can be used one-to-one with any of the alarm receivers or as part of a multipatient system.

ORDER CODE WPAS10 Price £61.53 (£73.83)



BED LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

CTM-3BKIT Bed leaving detection alarm with recordable voice alarm and pager

The majority of falls occur in the home when a person leaves their bed. With the CTM-3, a voice recording can be made for playback when a person rises from their bed, along with an alarm transmission to the included MPPL tone/vibrate carer alarm pager.

- Can help reduce falls and potential wandering with advanced warning of bed leaving
- Integral transmitter with operating range of 100M
- Under mattress pressure sensor 750mm x 450mm
- Instant (3 second) or 15 minute delayed alarm
- CTM-3 can switch to floor pressure mat mode
- CTM-3 Powered by 3 x AAA batteries or AC adaptor



Bed leaving kit ORDER CODE CTM3-BKIT PRICE £94.00 (£112.80)

NMDRX-CTMB Bed leaving detection alarm with wireless audible/flashing alarm

This kit includes the CTM-3 voice alarm transmitter, bed pressure mat CT-3 and the NMDRX wireless tabletop alarm. The NMDRX is a 5-Channel alarm, which enables the inclusion of additional sensors, such as: Door alarms, worn fall sensors (FA6) or PIR (TUMPIR21). A comprehensive home care system is easily configured.

- NMDRX suitable for deaf/hearing impaired carers
- 5 Independent flashing LED colours to assign to alarm channel



Complete kit ORDER CODE NMDRX-CTMB PRICE £102.00 (122.40)

NMDTXRXM-CT3B Bed leaving detection kit with channel selectable alarm transmitter

This kit includes the CT-3 bed pressure mat, NMDRX 5-Channel alarm receiver and the NMDTXM multifunction sensor transmitter. The transmitter channel codes are easily changed to allow additional bed monitors to be included in a system (up to 5), with the NMDRX indicating an active alarm with different coloured flashing channel LED's.

- Med Page
- NMDTXM powered by 2 x AAA batteries
- Transmission range up to 200M

Complete kit ORDER CODE NMDTXRXM-CT3B PRICE £96.31 (£115.57)

TUMAC1 Home carer falls prevention kit

TumbleCare®

The TUMAC1 Provides an effective solution for detecting and dealing with falls quickly in the home. During the day, the fall sensor can be worn with the included wrist strap or lanyard. The gravity change sensor responds rapidly to sudden changes in height resulting in an alarm transmission to the caregiver pager. For night-time monitoring the bed alarm pressure sensor is placed under the bed mattress. If the user exits the bed an alarm signal is transmitted to the caregiver pager.

Key Features

- An effective solution to detect and help prevent falls in the home or in care
- Includes fall detection sensor pendant and bed leaving alarm for round the clock monitoring
- Waterproof fall sensor with SOS call button ideal for showering
- Caregiver pager with tone and silent vibration for effective alarm notification



This kit includes: 1 x TXP11 fall pendant, 1 x TUMPAG31 caregiver pager, 1 x BTXMED04 bed alarm sensor monitor, 1 x MED04BP bed sensor pad, batteries included.

ODER CODE TUMAC1 PRICE £127.09 (£152.51)

TUMWPASKIT Cordless bed sensor pressure mat with caregiver pager

 New central and edge detect pressure sensor pad with concealed transmitter

 Cargiver pager with tone & vibration alarm notification (2 x AA batteries)

 Raises an alert to carer when patient sits on or falls off edge of bed

 Pager can monitor 3 x bed sensors (see WPAS10 page 29)

Ready to use straight out of the box



MADE SIMPLE

ORDER CODE TUMWPASKIT PRICE £89.73 (£107.68)

BED LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

CTMV-CT3 Bed leaving detection alarm with recordable voice prompt

This kit is suited to people who may be confused when waking during the night. Upon rising from the bed a personalised message is played back as a warning or instruction. The CTMV is a receiver with 3 voice alarm channels that can be triggered by remote transmitters. A door alarm, for example, could remind a person a door is open before getting into bed.



- Under mattress pressure sensor 750mm x 450mm
- 3-Channel voice recordings of up to 15 seconds each
- Voice alarm CTMV powered by 3 x AAA batteries (included)

Complete kit ORDER CODE CTMV-CT3 PRICE £56.83 (£68.20)

BMDV11-BKIT Dual alarm bed monitor with recordable voice prompt and clothing clip pull-away alarm

This multi-function bed/chair exit alarm monitor from Rondish, offers state-of-the-art fall prevention technology for home and professional care use. The dual alarm trigger function (sensor pad or clothing clip) provides options for bed,chair or W/C. Can be cable connected to Nurse Call system.

- Pull cord clip alarm use for bed,chair, W/C
- Bed under-mattress pressure sensor
- Recordable voice prompt "stay in bed" etc

Complete kit ORDER CODE BMDV11-BKIT PRICE £82.12 (£98.54)



BMS-03BKIT Bed leaving alarm under-mattress pressure sensor, audible alarm and nurse call connection output

Basic bed leaving alarm which produces an audible alarm when a patient rises from their bed.

- Under mattress pressure sensor
- Sensor controller with audible alarm
- Sensor pad fault/unplugged alarm
- Powered by 3 x AAA batteries
- Nurse Call connection output jack socket

Rondish, San American Company (Control of Control of Co

Controller with sensor ORDER CODE BMS-03BKIT PRICE £66.00 (£79.20)

BED LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

TumbleCare by Medpage bed exit detection with pager alarm



The bed sensor is positioned under a patient's bed mattress to detect their presence or exit. On exit the alarm control transmitter sends a radio alarm signal to the carer radio pager to produce an audible and vibrating alarm notification.

- Complete bed exit detection alarm system
- Radio pager notification to bed exit
- Reduce falls in domestic and professional care
- Hard wearing materials
- Rubber protection jacket for alarm transmitter

ORDER CODE TUMBSMPPLK PRICE £100.30 (120.36)

CBM-03PAGBKIT Bed leaving alarm kit with wire-free under-mattress pressure sensor and caller ID display pager

Totally cordless bed leaving alarm system. The alarm controller transmitter is concealed in the corner of the sensor. With no visible signs of the sensor under the mattress, tampering and damage by wire pulling is reduced. Kit includes PAG-11C caller ID pager and CBP-01 bed sensor with transmitter.

- Suitable for single or multiple bed monitoring
- Compatible with CMU-02, CMD-11C Nurse Call stations

Kit price ORDER CODE CBM-03PAGBKIT PRICE £131.30 (£157.56)

NMDTX-BKIT Under mattress bed leaving sensor for Care Homes

This kit is recommended for larger buildings, stone-built properties or buildings with known signalling difficulties. The bed pressure mat is 750mm x 450mm and is positioned under the bed mattress. When the patient rises from their bed a data encoded display message is transmitted to the carer pager. For example, 'Bed 1'.



The kit comprises of NMDTX2 long range POCSAG transmitter, 1 x CT-3 pressure mat and 1 x EM300 pager. An unlimited number of long range devices can be configured into a multiple care system. Single kit default message "Bed Alarm".

Single user kit ORDER CODE NMDTX-BKIT PRICE £145.00 (£174.00)

CHAIR LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

CTM-3CKIT Chair leaving detection alarm with recordable voice alarm and pager

Falls often occur in the home when a person rises from their chair. With the CTM-3, a voice recording can be made for playback when a person rises from their chair, along with an alarm transmission to the included MPPL tone/vibrate carer alarm pager.

- Can help reduce falls and potential wandering with advanced warning of chair leaving
- Integral transmitter with typical range of to 100M
- Under cushion pressure sensor 380mm x 300mm
- Powered by 3 x AAA batteries or optional AC adaptor
- Monitor holding bracket included

Complete kit ORDER CODE CTM-3CKIT PRICE £90.75 (£108.90)

NMDRX-CTMCK Chair leaving detection alarm with wireless audible/flashing alarm

This kit includes the CTM-3 recordable voice alarm transmitter, chair pressure mat CT-2 and the NMDRX wireless tabletop alarm. The NMDRX is a 5-Channel alarm, which enables the inclusion of additional sensors, such as: Door alarms, worn fall sensors (FA6) or PIR (TUMPIR21). A comprehensive home care system is easily configured (See page 12).

- NMDRX suitable for deaf/hearing impaired carers
- 5 Independent flashing LED colours to assign to alarm channel - 100M wireless range
- Voice reminder and alarm receiver alert on chair exit





Complete kit ORDER CODE NMDRX-CTMCK PRICE £99.99 (£119.99)

NMDRX-NMDTXCKITChair leaving detection kit with channel selectable alarm transmitter

This kit includes the CT-2 chair pressure mat, NMDRX 5-Channel alarm receiver and the NMDTXM multi-function sensor transmitter. The transmitter channel codes are easily changed to allow additional chair monitors to be included in a system (up to 5) with the NMDRX indicating an active alarm with different coloured flashing channel LED's.

- NMDTXM powered by 2 x AAA batteries
- Transmission range up to 100M
- Advanced warning fall prevention/wander detection
- Portable alarm receiver with audible/flashing LED's

Easylink UK

Final and the second of the sec

Complete kit ORDER CODE NMDRX-NMDTXCKIT PRICE £88.81 (£106.57)

CHAIR LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

CTMV-CT2 Chair leaving detection alarm with recordable voice prompt

This kit is suited to people who may be unsteady on their feet when rising from their chair. A personalised recordable message is played back as a warning or instruction.

- Under cushion pressure sensor 380mm x 300mm
- 3-Channel voice recordings of up to 15 seconds each
- Integral alarm receiver with Nurse Call connection output
- Voice alarm CTMV powered by 3 x AAA batteries (included)
- · Advanced warning of potential fall or wandering
- Single or looped message playback

Complete kit ORDER CODE CTMV-CT2 PRICE £56.83 (£68.20)

BMDV-11CKIT Dual alarm chair monitor with recordable voice prompt and clothing clip pull-away alarm

This multi-function chair exit alarm monitor from Rondish, offers state-of-the-art fall prevention technology for home and professional care use. The dual alarm trigger function monitor provides options for under chair cushion (pad) or pull-off clothing clip, also useful for W/C. Can be cable connected to Nurse Call.



- Pull cord alarm attaches to clothing
- Under chair cushion pressure sensor
- Recordable voice prompt "stay in your chair" etc.

Complete kit ORDER CODE BMDV-11CKIT PRICE £78.40 (£94.08)

BMS-03CK Chair leaving alarm under-cushion pressure sensor, audible alarm and nurse call connection output

Basic chair leaving alarm which produces an audible alarm when a patient rises from their chair.

- Under-cushion pressure sensor
- Sensor controller with audible alarm
- Volume adjustment
- Sensor pad fault/unplugged alarm
- Powered by 3 x AAA batteries
- Nurse Call connection output jack socket

Rondish, grave assert across

Controller with sensor ORDER CODE BMS-03CK PRICE £62.00 (£74.40)

CHAIR LEAVING ALARMS FOR SINGLE OR MULTIPLE PATIENTS

TumbleCare by Medpage chair occupancy detection alarm system

The chair sensor is positioned under a patient's chair cushion to detect their presence or exit. On exit the alarm control transmitter sends a radio alarm signal to the carer alarm station to produce an audible and flashing light alarm notification.

- Complete chair exit detection alarm system
- Wireless carer alarm notification to chair exit
- Reduce falls in domestic and professional care
- Hard wearing materials
- Quality chlorine proof chair sensor
- Suitable for single or multiple patient chair monitoring (with additional monitors)



5-Channel alarm station

ORDER CODE TUMCSDRXTK PRICE £100.30 (120.36)

CBM-03PAGCKIT Chair leaving alarm kit with wire-free under-cushion pressure sensor and caller ID display pager

Rondish

Totally cordless chair leaving alarm system. The alarm controller transmitter is concealed in the corner of the sensor. With no visible signs of the sensor under the cushion, tampering and damage by wire pulling is reduced. Kit includes PAG-11C caller ID pager and CCP-01 chair sensor with transmitter.

- Expandable from a single to multiple patient system
- · Alarm signal transmitted to pager upon rise from chair
- Compatible with CMU-02 & CMD-11C nurse call stations

Kit price ORDER CODE CBM-03PAGCKIT PRICE £101.10 (£121.32)

NMDTX-CKIT Under-cushion chair leaving alarm for Care Homes

This kit is recommended for larger buildings, stone-built properties or buildings with known signalling difficulties. The chair pressure mat is 380mm x 300mm and is positioned under the chair cushion. When the patient rises from their chair a data encoded display message is transmitted to the carer pager. For example, 'Chair Alarm 1'.

The kit comprises of NMDTX long range POCSAG transmitter, 1 x CT-2 pressure mat and 1 x EM300 pager. An unlimited number of long range devices can be configured into a multiple care system. Single kit default message "Chair Alarm".

Single user kit ORDER CODE NMDTX-CKIT PRICE £139.00 (£166.80)



MEDAD10P Waterproof wireless SOS pendant with telephone auto dialler

The MEDAD10 automated dialler uses the mobile telephone network to automatically call up to 5 trusted stored emergency contacts. Pendant activation will trigger the dialler into making telephone calls to nominated phone numbers allowing 2-way voice calling as well as sending an emergency SMS text message

SOS Both Ny man

Features

- Perfect for the elderly and vulnerable people living alone
- Waterproof SOS call button shower and bath safe
- Button transmitter range 100M
- Mobile network GSM SOS emergency telephone dialler
- Completely automated emergency contact dialling
- Setup and alert receiving APP for iOS and Android
- No monitoring or subscription charges
- Simple system setup by telephone keypad (including land line telephone) or SMS Text messages

Technical

- Integral radio receiver to pair with up to 30 transmitters
- Label transmitters for notification ID (first 5)
- SOS one-touch emergency dialling
- 2 x Wired trigger inputs N/O N/C configurable
- Audio listen-in function from remote phone
- Comprehensive system status monitoring
- Power: DC 12V power adaptor (included)
- GSM bands: 850/ 900/ 1800/ 1900MHz
- Unlocked to any mobile network

ORDER CODE MEDAD10P PRICE £74.99 (89.99)

1PMobile GSM SIM card - the best and cheapest way to run GSM based products

You can purchase your 1P Mobile SIM card through Medpage with your GSM device, tracker or security dialler to receive 1p a minute phone calling, 1p a text, 1p a MB data. Join 1pMobile today with a £30 top-up so you are all set for a whole year, or £10 for 4-months. Our triple SIM card is free, and it will work in any unlocked 2G, 3G or 4G phone or GSM device.



Best mobile network coverage (UK) Stay connected in known poor signal areas Ideal for GPS trackers and security diallers

SEE THIS PRODUCT ON OUR WEBSITE